

Premium quality and comfort...

at a budget- friendly price.

**The perfect treatment table for a
small space and a limited budget.**

aXionXelect requires only 15 square feet of space, thanks to its vertical-raising scissor mechanism.

Its adjustable sections ensure the table remains stable throughout and simple to operate. The aXionXelect treatment table is a cost-effective option for a practice with little available space.

Zimmer
MedizinSystems



aXionXelect

**In stock and ready to
ship immediately!**




- Electrical height adjustment from 19" to 41.5"
- Section adjustment with gas springs
- Ergonomic handles for operation with one hand
- Scissors movement saves space - the treatment table moves up and down linearly
- Padded sections with soft edges
- Optimal patient comfort, even when lying face-down, is ensured by the face hole with triangular insert.

**Quality you are accustomed to with a
Zimmer MedizinSystems product.**


axionXelect



Optional:
1550800 Paper roll holder for 50 cm rolls, mounted under foot-section




axionXelect




Our imitation leather upholstery is flame retardant and phthalate free –

A unique combination that makes this material suitable for medical use!



- Tabletop size (L x W): 81.5" x 26.4"
- Height adjustment: approximately 19" - 41.5"
- Lifting time (min.-max.): approximately 35 sec.
- Lifting capacity: 550 lbs
- Power supply: 115 volts
- Color code frame: Grey = RAL 7035
- Color code upholstery: Grey 213
- Foam density: 2.4 lb/ft³



Removable wheel set can be raised and lowered for easy moving - and locked for safety and stability when in use for treatment.

Zimmer MedizinSystems
3 Goodyear, Suite B
Irvine, CA 92618
(800) 327-3576 (office)
info@ZimmerUSA.com
www.ZimmerUSA.com
shop.ZimmerUSA.com

Zimmer
MedizinSystems