

**Versatile
positioning.**

**Professional
comfort.**

Z Ultimate



3-Section Treatment Table

**with adjustable arms
and Ergo-Wing segments**

The Z Ultimate 3-Section by Zimmer MedizinSystems is a premium electric treatment table designed for modern therapy clinics, combining flexible positioning, smooth height adjustment, and a stable rolling lift system to support efficient daily workflows and patient comfort.

**Quality you are accustomed to with a
Zimmer MedizinSystems product.**



Z Ultimate

- ✓ Adjustable arm sections with integrated Ergo-Wing support
- ✓ Thick, comfortable padding with soft edges
- ✓ Stable frame with integrated Autolock safety system
- ✓ Durable construction designed for high-use clinical environments



Key Features

- Electric height adjustment
- Gas-spring assisted section positioning
- Fold-away head section for cervical access
- Ergonomically positioned control handles
- Rolling lift system for easy mobility
- Adjustable arms and Ergo-Wing segments

Technical Data

Table Surface:	Approx. 80.3" × 26.4"
Height Range:	18.5" – 39.4"
Lift Time:	Approx. 18 sec
Capacity:	440 lbs
Power:	120 VAC, 60 Hz
Color Code Frame:	Anthracite (RAL 7016)

Zimmer MedizinSystems
3 Goodyear, Suite B
Irvine, CA 92618
(800) 327-3576 (office)
info@ZimmerUSA.com
www.ZimmerUSA.com
shop.ZimmerUSA.com

Zimmer MedizinSystems